# Caframo Real Torque Digital Stirrers

## Power and Control

ENGLISH

Instruction Manual  
pages 2-17

FRANÇAIS

Mode d=emploi  
pages 18-33

ESPAÑOL

Manual de instrucciones  
páginas 34-49

DEUTSCH

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RoHS  
WEEE

 **Caframo**<sup>®</sup>

# CONTENTS

Safety Instructions .....	3
General Instructions .....	4
Stirrer Approvals .....	4
Packing List .....	5
Description .....	6
Technical Data .....	7
BDC Torque Capabilities .....	8
BDC Specifications .....	9
BDC Keypad Membrane .....	9
Installation in the Laboratory .....	10
Environmental Operating Conditions .....	10
General BDC Operating Instructions .....	11
Routine Cleaning, Service, Transport & Storage .....	13
Warranty and Product Liability .....	14
Caframo Accessories .....	15

**On-Line Warranty Registration**

**[www.caframo.com/warranty.htm](http://www.caframo.com/warranty.htm)**

# Safety Instructions

-  1. All operators must be familiar with the stirrer and should read this entire manual.
-  2. The stirrer must be securely fixed to a stable support. Mount to a Caframo Stand (Part Number A110) with a Caframo heavy duty clamp (Part Number A120). If other stands or clamps are used, the unit must be secure and can not fall if tipped up to ten degrees from the vertical.
-  3. The stirrer must be supplied with the rated voltage. See serial plate.
-  4. CAUTION: This is not an explosion proof stirrer. Do not use with highly flammable or explosive materials.
-  5. Spinning paddles or impellers can cause severe personal injuries. Operators must use extreme care and good judgement when mixing at any speed. The BDC stirrers have higher torque capabilities than conventional stirrers.
-  6. All mixing paddles and impellers must be in good condition with straight shafts. If the stirrer vibrates at high speeds check the paddle shaft for damage and repair or replace it.
-  7. Extreme care must be taken when mixing chemicals so that no chemicals are splashed outside the mixing vessel. Care must be taken when changing to faster mixing speeds. On power up, the unit will display its set speed and will climb to that speed when the pause button is pushed. Always start at lowest speed if unsure of maximum safe speed.
-  8. Ensure that the mixing impeller does not contact the containing vessel.
-  9. Do not operate while standing in water. Keep the unit dry and do not immerse any part, except the mixing paddle, into any liquids. Protect from splashing.
-  10. Ensure that no loose clothing, jewelry, or hair can become entangled in any rotating parts. Fast spinning chuck can cause injury to operator. Use chuck guard provided when stirrer is in use.
-  11. Power can be interrupted to the stirrer by pressing the power button or by disconnecting the mains cord. If rotating, this will cause the stirrer to stop and will disconnect power to the stirrer's internal drive circuit.
-  12. Shaft rotation can also be stopped by pressing the pause button. This does not disconnect power to the stirrers internal drive circuit.
-  13. Wear safety goggles and suitable clothing when operating the stirrer.
-  14. Repairs must be carried out only by technicians authorized by Caframo.

# General Instructions

1. Read all safety instructions.
2. Unpack the instrument carefully. Ensure the instrument is not damaged and verify the contents against the packing list. If the instrument is damaged or the contents do not correspond with the packing list, notify your supplier immediately.
3. Read the manual carefully! Make sure that all users read this manual thoroughly before operating the BDC stirrer.

4. This instruction manual should always be kept readily available for reference.
5. For further questions or inquiries, please contact:

## **Caframo Ltd.**

R.R. # 2, Warton, Ontario, Canada N0H 2T0  
 Tel: (519) 534-1080 Fax: (519) 534-1088

### Toll free in USA and Canada:

Tel: (800) 567-3556 Fax: (800) 209-6786

Web Site: [www.caframo.com](http://www.caframo.com)

E-mail: [caframo@caframo.com](mailto:caframo@caframo.com)

## Stirrer Approvals

Models	Standards
BDC1850 BDC3030 BDC6015	 CAN / CSA C22.2 (1010-1) <b>IEC1010</b>
BDC1850-220 BDC3030-220 BDC6015-220	 <b>CE, IEC1010, EMC EN50081-1,            EN50082-1 (EN61000-4-2,            EN61000-4-3, EN61000-4-4,            EN61000-4-5)            RoHS            WEEE</b>

CSA Canadian Standards Association

IEC Low Voltage Act

EMC Electromagnetic Compatibility Generic Immunity Standard

RoHS Restriction of Hazardous Substances (EU)

WEEE Waste Electrical and Electronic Equipment (EU)

# Packing List

Compare the contents of the shipment to the items on the following list to ensure all parts are received with the unit. Do not discard the container and packing materials until all parts are accounted for. It is recommended to keep the original packing for transportation and storage.

Description	Order Number	Quantity
Ultra Torque Stirrer, 115 volt	BDC1850	
or Universal Stirrer, 115 volt	BDC3030	
or Ultra Speed Stirrer, 115 volt	BDC6015	1
or Ultra Torque Stirrer, 220 volt	BDC1850-220	
or Universal Stirrer, 220 volt	BDC3030-220	
or Ultra Speed Stirrer, 220 volt	BDC6015-220	
Calibration Certificate	BDC00PK013	1
Chuck with key	13C-SET	1
Chuck Guard	CG-1	1
Instruction Manual	BDCMAN rev05	1

## Description

All Caframo stirrers are specifically designed for professional use in the laboratory and for use in industry to stir numerous liquids with various properties and applications. The BDC series of stirrers are available in six different models.

BDC1850  
BDC1850-220            rotates at 12 - 1800 rpm with a maximum torque of 568 Ncm (50 in-lbs)

BDC3030  
BDC3030-220           rotates at 20 - 3000 rpm with a maximum torque of 341 Ncm (30 in-lbs)

BDC6015  
BDC6015-220           rotates at 40 - 6000 rpm with a maximum torque of 170 Ncm (15 in-lbs)

\* The “220” models operate on 220-240 volt mains voltage.

The BDC stirrers use a digitally controlled high efficiency brushless DC motor which yields high torque over a very wide speed range. The speed is kept constant with continuous feedback to the motor. The BDC stirrers are designed for reliability and durability. There are no wear parts that require user service during the normal life of the instrument. There are no user serviceable parts.

The BDC stirrers incorporate smart logic for automatic transmission control. This logic determines the appropriate speed range at which the stirrer should operate. See graphs on page 8. The automatic selection can be overridden by pressing the stirrer range button.

With higher solution viscosities, more torque is required to stir at the same speed. When the torque demands are outside the instrument’s range, the control logic slows down the stirrer until the torque comes within the instrument’s range. If the torque demand is too great at the stirrer’s lowest possible speed, the unit will automatically shut down.

The BDC stirrers are designed for safety. The unit comes on in “pause” mode once power is activated. After a power interruption the stirrer remains off. The stirrer’s speed will gradually ramp up to the desired set point after the pause button is pressed.

All BDC stirrers have integrated electronic overload protection which ensures no damage can occur to the electronics, motor or transmission when an overload situation occurs, such as too great a torque demand or unit overheating.

A chuck guard is supplied to place over the chuck. It protects the user from the spinning chuck and also helps to protect the chuck from potential splashes.

# Technical Data

# BDC1850 / 3030 / 6015

Voltage	120 volts @ 60 Hz	230 volts @ 50/60 Hz
Installation Category (in accordance with IEC664)	II	II
Maximum Current Draw	3 Amps	3 Amps
Maximum Output Power	1/5 hp, 150 watts	
Speed Readout	LCD 4 digit display	
Speed Readout Accuracy	+/- 1 rpm or +/- 1% of reading	
Torque Readout Accuracy	+/- 12 N-cm (1 in-lb) or +/- 5% of reading	
Maximum Operating Noise	70 dB(A) at 1 meter	
Motor Type	Brushless DC	
Drive	2 stage transmission	
Hollow drive shaft	Nitrated P .30 with 10 mm (3/8") diameter bore for stirring impellers 13 mm (0.51") outside diameter for drive chuck	
Chuck Guard	Length 67 mm, (2.64") Diameter 54 mm, (2.13") Molded of clear chemical resistant polyurethane.	
Seals	Butyl Rubber IP 42	
<b>Holding chuck</b>	<b>3 prong stainless steel. Maximum opening 9.5 mm (3/8") diameter</b>	

Protection measures:		0
Class	Protection device	Recommended use
0	No additional protection	Used only on equipment where failure of the controls does not represent a hazard.

**The operating conditions should be monitored at regular intervals.**

## BDC Torque Capabilities

Stirrer	Low Speed			Cont. Torque* Ncm	High Speed			Cont. Torque* Ncm
	Min Speed rpm	Max Speed rpm	Peak Torque Ncm		Min Speed rpm	Max Speed rpm	Peak Torque Ncm	
<b>BDC1850</b>	12	350	568	320	124	1800	114	75
<b>BDC3030</b>	20	584	339	190	206	3000	68	43
<b>BDC6015</b>	40	1168	170	85	412	6000	34	21

\*Continuous torque at 50% of rpm range and 25 degrees Celcius.

Torque Conversion Table		50 in-lb = 800 in-oz = 568 N-cm
In-lb to In-oz	Multiply by 16	30 in-lb = 400 in-oz = 339 N-cm
In-lb to N-cm	Multiply by 11.36	15 in-lb = 240 in-oz = 170 N-cm
N-cm to In-oz	Multiply by 1.416	

Fig 1


Fig 2



Fig 3

**BDC Specifications**


Weight: 5 kg (11 lbs)  
 Mounting Rod: 12.7 mm (1/2") diameter  
 140 mm (5-1/4") long

Display Accuracy

Speed: ± 1 rpm or ± 1% of reading  
 Torque: ± 12 N-cm or ± 5% of reading

**Keypad Membrane**

Fig 4

**Function Buttons**

- SPEED
- DISPLAY
- POWER
- PAUSE
- ± (RPM)

**Indicator LEDs**

- 1 High Speed
- 2 Low Speed
- 3 RPM
- 4 Torque (N-cm)
- 5 Overload
- 6 Power
- 7 Pause

# Installation in the Laboratory

1. The stirrer must be securely fixed to a stable support - e.g. clamped to an original Caframo stand that is on a stable, level surface. See safety instructions on page 3.

## CAUTION:

This is not an explosion-proof stirrer. Do not use with highly flammable or explosive materials.

## WARNING:

Be sure that the unit's voltage is the same as the supplied voltage. Information on the unit's voltage and current are located on the serial plate.

This unit must always be connected to a grounded electrical outlet.

2. Chuck removal and attachment:

Loosen the set screw with the allen key found on the chuck key. Remove chuck. To attach align the set screw with the pilot hole on the shaft. Tighten the set screw with the allen key found on the chuck key. Make sure that the set screw is firmly tightened. Open the chuck fully to prepare for impeller attachment.

3. Impeller attachment:

Insert the impeller shaft (see accessories) into the open chuck and adjust the impeller to a maximum length of 25 cm (10 inches) [of a 40 cm (16 inch) shaft] to reduce vibration. The hollow shaft allows the impeller shaft to come through the top of the stirrer if necessary. Tighten the chuck's outer ring by hand and then secure with the chuck key. Install the chuck guard over the chuck.

# Environmental Operating Conditions

All BDC stirrers must operate in the following conditions:

## Indoors

1. Altitudes below 2000 meters.
2. Temperatures from 5 °C to 40 °C
3. Maximum relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40 °C.
4. Mains supply as indicated under the technical specifications.
5. Pollution degree 2 in accordance with IEC 664

“Normally only non conductive pollution occurs. Occasionally, however, temporary conductivity caused by condensation must be expected”

# General BDC Operating Instructions

-  1. Make sure that the stirrer is securely mounted as specified in the installation instructions. Also make sure that the mixing paddle or impeller is securely attached in the desired position.
-  2. Connect the power cord to a grounded supply.
-  3. The stirrer is now ready for use to mix liquids and liquid/solid solutions. The operator, ensuring that all safety instructions are followed must determine appropriate impeller and container specifications.
-  4. The stirrer can operate continuously if the torque demand falls below the continuous torque operation line (see graphs on page 8). Should the ambient temperature be above 25 °C, then the 40 °C continuous operation line must be used.
-  5. If the stirrer is used for an extended period of time above the appropriate continuous temperature operation line (see graphs on page 8), then it may shut down. The stirrer may be started after it has cooled down.

## Button Operating Functions

- POWER - when power button is pressed for the first time stirrer is activated in “pause” mode. The set speed is displayed. RPM LED will be flashing indicating this is a set speed.
- When stirrer is on the POWER LED is on.
- When the POWER button is pressed a second time, the stirrer is turned off. The last set speed is remembered.
- PAUSE - stops the stirrer when operating.
- PAUSE returns the stirrer to the previous condition and previous display or starts the stirrer if stopped.
- PAUSE LED is on when stirrer is paused and set speed is displayed.
- DISPLAY - When running, toggles display between speed (RPM) and torque (N-cm). The LEDs indicate which parameter is displayed.
- When paused, the display will always show set speed.
- SPEED - determines speed range (high speed/low speed) and overrides automatic choice when allowable. The LEDs indicate which speed range is operating.
- + / - BUTTONS - Press buttons once to increase / decrease set speed by one rpm. Press and hold to increase / decrease speed at a greater rate.
- When unit is in pause, user can change the set speed by pressing the buttons.
- RPM LED will flash when displaying desired set speed and will be constant when displaying actual speed.

- As unit is auto ramping to a faster speed the user can stop the increase in speed by pressing the down arrow button. This changes the set speed to the present actual value.
- As unit is auto ramping to a slower speed the user can stop the decrease in speed by pressing the up arrow button. This changes the set speed to the present actual value.

### Indicator Lights

- HIGH SPEED / LOW SPEED LEDs - These indicate the stirrer's set range.
- OVERLOAD LED - This indicator will flash when an overload condition is approached. This indicator will remain on when an overload condition is reached and the stirrer is shut down.

Overload conditions will occur when:

- The torque demand is outside the unit's capabilities (solution is too viscous).
- The unit overheats due to extended operation at very high torque and/or speed demands.

When the overload condition is due to torque demand, the stirrer will automatically reduce speed to reduce the load. If reduction in speed is not sufficient to keep the load in the stirrer's capable range, then the stirrer will shut down.

When the overload condition is due to overheating, the stirrer will not automatically reduce speed. If the overload light is flashing and the speed has not reduced from the set point, the load should be reduced by the operator.

In either case, if the stirrer completely shuts down due to an overload condition, it can be reset by powering off the stirrer. If the stirrer has overheated then it must cool down before it can be reset.

# Routine Cleaning, Service, Transport & Storage

## Cleaning

Using a mild soap and water solution, wipe the housing and keypad of the stirrer with a damp cloth.

### Note:

Do not use chlorine bleach, chlorine-based cleanser, abrasives, ammonia, steel wool or scouring pads with metal content or similar harsh solvents or abrasives. These will damage the surface of the instrument.

## Service

For qualified repair technicians, disconnect power before servicing.

The BDC1850, BDC3030, BDC6015, are all fused with a “one-use” internal 3 amp fuse.

The BDC1850-220, BDC3030-220, BDC6015-220 are all fused with a “one-use” internal 3 amp fuse.

Do not attempt to service or repair a Caframo BDC stirrer. Contact Caframo for return authorization and instructions. See page 4.

## Transport and Storage

1. Place the unit and its parts into the original packing or a container with the necessary protection to prevent breakage during transport. Seal the original packing or container with packing tape.
2. Store the packed unit in a dry place.

### CAUTION:

When transporting this unit ensure it is not dropped or subject to unnecessary jolts.

Failure to clean, service, and handle the stirrer as outlined can cause the unit to fail and impair the safety of the user.

# Warranty and Product Liability

## Warranty

Caframo warrants the herein described product of its manufacture for three years from receipt of the equipment by the purchaser, against defects in material and workmanship. Warranty registration is available on-line ([www.caframo.com/warranty.htm](http://www.caframo.com/warranty.htm)) or by completing the supplied warranty registration card. Mail or fax back as indicated.

This limited warranty covers parts and labor. In the event of a warranty claim, contact Caframo. If the cause is determined to be a manufacturing fault, Caframo will repair or replace all defective parts to restore the unit to operation.






This warranty may only be altered by a specifically published amendment. No individual has authorization to alter the provisions of this warranty policy or its amendments. Expendable items are not covered by this warranty. Damages due to corrosion, accidental breakage or submersion are not covered.

## Product Liability






Under no circumstances shall Caframo be liable for indirect, consequential or special damages of any kind.

Caframo reserves the right to make technical changes without prior notice.

# Caframo Accessories

	Features	Caframo Part Number
	<p>Safety Stand</p> <ul style="list-style-type: none"> <li>base cast zinc alloy with chemical-resistant epoxy coating</li> <li>includes 710 mm (28") rod, stand base and allen key</li> <li>tube-stainless steel 304, Ø 25.4 mm, (1")</li> <li>5.0 kg, (11 lbs)</li> </ul>	<p>A110</p> <p>ROD-28 28" rod only ROD-48 48" rod only A110BASE base only</p>
	<p>Compact Stand</p> <p>to be used with CompactDigital™ Stirrer</p> <ul style="list-style-type: none"> <li>base cast zinc alloy with chemical-resistant epoxy coating</li> <li>includes 710 mm (28" ) rod, stand base and allen key</li> <li>rod - solid stainless steel 304, Ø 16 mm (5/8")</li> <li>5.0 kg, (11 lbs)</li> </ul>	<p>A210</p>
	<p>Fastening Clamp</p> <ul style="list-style-type: none"> <li>Used with A110 or A210 above</li> <li>Accepts 15 mm Ø through 30 mm Ø rods (0.6" - 1.2")</li> <li>Supports stirrer on the rod of the stand</li> <li>Cast zinc alloy with chemical-resistant epoxy coating</li> <li>Stainless steel 304 handles</li> <li>640 g, (12 lbs)</li> <li>carton quantity of 4</li> </ul>	<p>A120</p>
	<p>Crossed Blade Impeller</p> <ul style="list-style-type: none"> <li>welded to shaft</li> <li>stainless steel 316L</li> <li>weighs 128 grams</li> <li>very high shear with high vortexing characteristics</li> <li>will vortex 5 gallons of water at 900 rpm</li> </ul>	<p>A130</p>
	<p>Straight Blade Impeller</p> <ul style="list-style-type: none"> <li>welded to shaft</li> <li>stainless steel 316L</li> <li>weighs 124 grams</li> <li>high shear with moderate vortexing characteristics</li> <li>will vortex 5 gallons of water at 1200 rpm</li> </ul>	<p>A131</p> <p>A231 Compact version of A131. Blade diameter is 31 mm (1.25")</p>


# Caframo Accessories

Impeller	Features	Caframo Part Number
	<p>Pivoting Blade Impeller</p> <ul style="list-style-type: none"> <li>pinned to shaft</li> <li>stainless steel 316L</li> <li>weighs 126 grams</li> <li>moderate shear with moderate vortexing characteristics</li> <li>will vortex 5 gallons of water at 800 rpm</li> </ul>	A140
	<p>Collapsible Blade Impeller</p> <ul style="list-style-type: none"> <li>pinned to shaft</li> <li>stainless steel 316L</li> <li>weighs 127 grams</li> <li>high shear radial flow for mixing in carboys or narrow neck flasks</li> <li>collapsed blades enter any opening over 10 mm diameter</li> <li>will vortex 5 gallons of water at 1000 rpm</li> </ul>	A141
	<p>Square Blade Paddle</p> <ul style="list-style-type: none"> <li>welded to shaft</li> <li>stainless steel 316L</li> <li>weighs 180 grams</li> <li>for large volumes or high viscosities; axial flow with maximum vortexing</li> <li>will vortex 5 gallons of water at 450 rpm</li> </ul>	<p>A150</p> <p>A250</p> <p>Compact version of A150. Blade diameter is 38 mm (1.5")</p>
	<p>Pitched Blade Propeller</p> <ul style="list-style-type: none"> <li>set screws fasten propeller to shaft</li> <li>stainless steel 316L</li> <li>weighs 181 grams</li> <li>low shear, high top-to-bottom pumping capacity with high gas injection</li> <li>will vortex 5 gallons of water at 3600 rpm</li> <li>propellers sold separately; can add multiple propellers to one shaft 8mm (5/16") Ø</li> </ul>	<p>A160</p> <p>A161 (propeller only)</p> <p>A162 (rod only)</p>
	<p>Baffle</p> <ul style="list-style-type: none"> <li>to fit a 4L beaker (custom sizes available)</li> <li>stainless steel 316L</li> <li>teflon bumper</li> </ul>	A170
	<p>Chuck Guard</p> <ul style="list-style-type: none"> <li>included with each stirrer</li> <li>clear chemical resistant polyurethane</li> <li>magnetic, hinged closure</li> </ul>	CG-1
	<p>Chuck</p> <ul style="list-style-type: none"> <li>3 prong stainless steel.</li> <li>maximum opening 9.5 mm (3/8") diameter</li> <li>includes key</li> </ul>	13C-SET

U Shaped (Anchor) Paddles

- Made of stainless steel 316L
- Blade, hub, and shaft welded together

Part #	Weight	A	B	C	D
U022	188 g (6.6 oz)	51 mm, (2")	51 mm, (2")	11 mm, (0.5")	11 mm, (0.5")
U044	280 g (9.9 oz)	101 mm, (4")	101 mm, (4")	25 mm, (1")	25 mm, (1")
U055	300 g (10.6 oz)	127 mm, (5")	127 mm, (5")	25 mm, (1")	25 mm, (1")
U510	716 g (25.3 oz)	127 mm, (5")	254 mm, (10")	50 mm, (2")	50 mm, (2")

\*U022 406mm (16"), Ø 8mm (5/16")

A600-SET

High Speed Collet and Shaft Set for high speed mixing (1800-6000 rpm)

Includes: dispersion blade, pitched blade propeller and a straight blade.

Total length of collet and shaft together: 26.4 cm (10 3/8")

Precision 5/8" diameter shaft of 316L stainless steel

Fits impellers with a 8 mm (5/16") diameter bore

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